



Apply water based adhesive like a pen with the **SOLO GLUE RITER™** from Gluefast!



Can replace double-sided tape!

The **SOLO GLUE RITER**

Press down on the item and adhesive flows!
Remove SOLO from the item and the flow stops!

Portable, versatile, and efficient.

The new **SOLO GLUE RITER** features a pen-like applicator attached to a simple squeeze bottle. Adhesive flows in 1/2" line widths.

- ✓ Designed for matboard gluing with many other uses!
- ✓ Many replace double sided tape with a 90% savings!
- ✓ No daily cleanup! When the **SOLO GLUE RITER** is not used the adhesive film dries but can be easily removed by peeling it away from the nozzle.
- ✓ The handy stand keeps the **SOLO GLUE RITER** ready for use with water based resin adhesives only! Dextrin glues are not recommended.
- ✓ Gluefast's **MATHESIVE 101** is fast setting and allows you to reposition the material for a few seconds. The **SOLO GLUE RITER** includes this adhesive with the bottle!

MATHESIVE 101 can be purchased in quarts, gallons, 5 gallon pails and 55 gallon drums. It is pH neutral (acid free). For more information request the Gluefast Service Bulletin for **MATHESIVE 101**.

THE P.A.T. JR. HAS A GLUE RITER WITH 1/8" OR 1/2" NOZZLES AND A PRESSURE TANK FOR LARGE VOLUME USERS. FOR SPECIFICATIONS, REQUEST LIT SHEET PAT-02.





Photo at Right: SOLO with stand



PART NUMBER	DESCRIPTION	NET WEIGHT	GLUE CAPACITY	BENCH SPACE
GRB/10-101	SOLO GLUE RITER with 10 ounce applicator & MATHESIVE 101	12 oz / 340 grams	10 ounces / 283 ml	2" Diameter